Patient Digital Image Delivery, Access, Viewing and Archiving: Recommended minimum standards for **DIAGNOSTIC PURPOSES** (PRIMARY or SECONDARY)

(© RACS DDI Options - www.surgeons.org/dirwp)

2. onto LCD monitor screen with or

without magnification

and size calibration

compliant format.



The (blue) boxes below represent current potential solution to address each process step.

To be to a standard suitable for diagnostic purposes, every green box must have a blue box solution. Any green box that does not have a blue box solution, appropriate to the particular clinical situation, potentially renders the images provided as not suitable for diagnostic use.